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1. Food and Drug Administration, 2000, Code of Federal Regulation Chapter 21.
2. Official Journal of the European Union, 2004, Regulation 1935/2004/EC.
3. Official Journal of the European Union, 2006, Regulation 2023/2006/EC.
4. Official Journal of the European Union, 2011, Regulation 10/2011/EU.
5. Official Journal of the European Union, 2006, Regulation 1907/2006/EC.
6. Official Journal of the European Union, 2009, Regulation 2009/49/EC.
7. Official Journal of the European Union, 1987, Directive 87/357/EEC.
8. Dzumana Z, Zachariasova M, Veprikovaa Z, Godulab M, Hajslova J. Multi-Analyte High Performance Liquid Chromatography Coupled to High Resolution Tandem Mass Spectrometry Method for Control of Pesticide Residues, Mycotoxins, and Pyrrolizidine Alkaloids. *Anal.Chim.Acta.* 2015; 863:29–40.
9. Pérez-Ortega P, Lara-Ortega FJ, García-Reyes JF, Gilbert-López B, Trojanowicz M, Molina-Díaz, A. A Feasibility Study of UHPLC-HRMS Accurate-Mass Screening Methods for Multiclass Testing of Organic Contaminants in Food, 2016, *Talanta* 160, 704–712.
10. Pérez-Ortega P, Lara-Ortega FJ, Gilbert-López B, Moreno-González D, García-Reyes JF, Molina-Día A. Screening of Over 600 Pesticides, Veterinary Drugs, Food-Packaging Contaminants, Mycotoxins, and Other Chemicals in Food by Ultra-High Performance Liquid Chromatography Quadrupole Time-of-Flight Mass Spectrometry (UHPLC-Q-ToF-MS). *Food Anal.Methods.* 2017; 10:1216–1244.
11. Romero-González R. Food safety: How Analytical Chemists Ensure It. *Anal.Methods.* 2015; 7:7193–7201.
12. Coscollà C, León N, Pastor A, Yusà V. Combined Target and Post-run Target Strategy for a Comprehensive Analysis of Pesticides in Ambient Air Using Liquid Chromatography-Orbitrap High

Resolution Mass Spectrometry.*J. Chromatogr.A.* 2014; 1368:132–142.

13. Sjerps RMA, Vughs D, van Leerdam JA, ter Laak TL, van Wezel AP. Data-Driven Prioritization of Chemicals for Various Water Types Using Suspect Screening LC-HRMS. *Water Research.*2016; 93:254–264.
14. Pringle SD, Giles K, Wildgoose, J. An Investigation of the Mobility Separation of Some Peptide and Protein Ions Using a New Hybrid Quadrupole/Travelling Wave IMS/oa-ToF instrument.*International Journal of Mass Spectrometry.*2014, 26, 1-12.
15. Giles K, Pringle SD, Worthington KR, Little D, Wildgoose J, Bateman RH.Applications of a Travelling Wave-Based Radio-Frequency Only Stacked Ring Ion Guide.*Rapid Commun.Mass Spectrom.*2004, 18, 2401.
16. McCullagh M, Pereira C.AM, Yariwake JH.Use of Ion Mobility Mass Spectrometry to Enhance Cumulative Analytical Specificity and Separation to Profile 6 - C/8 - C glycosylflavone Critical Isomer Pairs and Known–Unknowns in Medicinal Plants.*Phytochemical Analysis.* 2019, 30(4), 1–13.
17. McCullagh M, Douce D, Van Hoeck E, Gosciny S. Exploring the Complexity of Steviol Glycosides Analysis Using Ion Mobility Mass Spectrometry.*Anal.Chem.*2018, 90, 4585–4595.
18. McCullagh M, Pereira C.AM, Yariwake JH.Use of Ion Mobility Mass Spectrometry to Enhance Cumulative Analytical Specificity and Separation to Profile 6 - C/8 - C glycosylflavone Critical Isomer Pairs and Known–Unknowns in Medicinal Plants.*Phytochemical Analysis.*2019, 30(4), 1–13.
19. Gosciny S, McCullagh M, Far F, De Pauw E, Eppe G. Towards the Use of Ion Mobility Mass Spectrometry Derived Collision Cross Section as a Screening Approach for Unambiguous Identification of Targeted Pesticides in Food.*Rapid Commun.Mass Spectrom.*2019, 33(S2), 1–15.
20. Nye LC, Williams, JP, Munjoma NC, Letertre, MPM, Coen M, Bouwmeester R, Nicholson, JK, Plumb RS, McCullagh M, Gethings L, Lai S, Langridge J, Vissers JPC, Wilson ID.A Comparison of Collision Cross Section Values Obtained via Ion Mobility Spectrometry Following Direct Infusion and an Evaluation of U(H)PLC-IM-MS for the Characterisation of Metabolites in Rat Urine.*J Chromatogr A* 2019 Sep 27; 1602:386–396.
21. McCullagh M, Goshawk J, Eatough D, Mortishire-Smith R, Pereira CAM, Yariwake JH, Vissers JPC.Profiling of the Known-Unknown Passiflora Variant Complement by Liquid Chromatography - Ion

Mobility - Mass Spectrometry. *Vissers.Talanta* 221, 2021, 121311.

22. Righetti L, Dreolin N, Celma A, McCullagh M, Barknowitz G, Sancho JV, Dall' Asta C. Travelling Wave Ion Mobility-Derived Collision Cross Section for Mycotoxins: Investigating Interlaboratory and Inter-Platform Reproducibility. *J. Agric.Food Chem.* 2020, 68, 39, 10937–10943.
23. Goshawk J, Barknowitz G, and McCullagh M. A Workflow for Automatic MS Library Creation from Time-of-Flight Full-Spectra Data Processed in UNIFI.Waters, White Paper [720006783EN < https://www.waters.com/waters/library.htm?lid=135047403>](https://www.waters.com/waters/library.htm?lid=135047403) .2020 Feb.
24. Yang Y, Hu C, Zhong H, Chen X, Chen R, Yam KL.Effects of Ultraviolet (UV) on Degradation of Irgafos 168 and Migration of Its Degradation Products from Polypropylene Films.*J Agric Food Chem.* 2016 Oct 19;64(41):7866–7873.
25. Hermabessiere L, Receveur J, Himber C, Mazurais D, Huvet A, Lagarde F, Lambert C, Paul-Pont I, Dehaut A, Jezequel R, Soudant P, Duflos G. An Irgafos® 168 story: When the Ubiquity of an Additive Prevents Studying Its Leaching From Plastics.*Sci Total Environ.*2020 Dec 20; 749:141651.
26. Jungmin L, Sorbitol, Rubus Fruit, and Misconception, *Food Chemistry*, Volume 166, 2015, Pages 616–622.
27. Cheng YY, Yu JZ, Minimizing Contamination from Plastic Labware in the Quantification of C₁₆ and C₁₈ Fatty Acids in Filter Samples of Atmospheric Particulate Matter and Their Utility in Apportioning Cooking Source Contribution to Urban PM_{2.5}. *Atmosphere*, October 2020 11(10):1120.
28. Nordby H E, Nagy S. Fatty Acid Profiles of Orange and Tangor Juice Sac Lipids.*Phytochemistry*, 1971, Volume 10, Issue 3, Pages 615–619.
29. Lamine M, Gargouri M, Rahali FZ and Mliki A. Authentication of Citrus Fruits Through a Comprehensive Fatty Acid Profiling and Health Lipid Indices: A Nutraceutical Perspective.*Journal of Food Measurement and Characterization* 2019 Vol.13 Issue 3 Pages 2211–2217.
30. Khan, FA, *et al.* "Comparative Evaluation of Physiochemical and Gc-MS Analysis of Sour Oranges and Sweet Oranges Peels Oil." *Life Science Journal* 2013;10(10s), 205–209.

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[ACQUITY UPLC I-Class PLUS系统 <https://www.waters.com/134613317>](https://www.waters.com/134613317)

[离子淌度质谱 <https://www.waters.com/waters/134656158>](https://www.waters.com/waters/134656158)

[SYNAPT XS高分辨率质谱仪 <https://www.waters.com/waters/nav.htm?cid=135020928>](https://www.waters.com/waters/nav.htm?cid=135020928)

[UNIFI科学信息系统 <https://www.waters.com/134801648>](https://www.waters.com/134801648)

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© 2022 Waters Corporation. All Rights Reserved.

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