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2. Chen Y, Zhu Z, Qiu H, Lu S. Occurrence of Perchlorate, Chlorate, and Bromate in Drinking Water in Shenzhen and Related Human Exposure Risks. *Environmental Advances*. 2022. 2666–7657.
3. Constantinou.P, Louca-Christodoulou D, Agapiou A. LC-ESI-MS/MS Determination of Oxyhalides (Chlorate, Perchlorate, and Bromate) in Food and Water Samples, and Chlorate on Household Water Treatment Devices along with Perchlorate in Plants. *Chemosphere* 2019.235: 757–766.
4. Panseri S, Noble M, Arioli F, Biolatti C, Pavlovic R. Occurrence of Perchlorate, Chlorate, and Polar Herbicides in Different Baby Food Commodities. *Food Chemistry* 330. 2020.
5. Zhang W, Lu Q, Li Y, Hua Y, Zheng R. Occurrence and Exposure Assessment of Chlorate and Perchlorate in Food and Drinking Water from Fujian China.*Food Control* 153.2023.
6. Li M, Xiao M, Xiao Q, Chen Y, Guo Y, Son J, Li R, Li C, Zhu Z, Qiu H, Liu X, Lu S. Perchlorate and Chlorate in Breast Milk, Infant Formulas, Baby Supplementary Food and The Implications for Infant Exposure. *Environmental International* 158.2022.
7. De-Alwis J, Ross E, Adams S. 使用沃特世阴离子型极性农药分析色谱柱和UPLC-MS/MS测定婴儿奶粉中的氯酸盐和高氯酸盐.沃特世应用纪要.720006686ZH.2020.
8. QuPPE Methods for Food of Plant and Animal Origin.https://www.eurl-pesticides.eu/docs/public/tmpl_article.asp?CntID=887&LabID=200&Lang=EN <https://www.eurl-pesticides.eu/docs/public/tmpl_article.asp?CntID=887&LabID=200&Lang=EN>
9. Document No.SANTE/11312/2021.Guidance Document on Analytical Quality Control and Method Validation Procedures for Pesticide Residues Analysis in Food and Feed.2021.

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waters_connect <<https://www.waters.com/waters/nav.htm?cid=135040165>>

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© 2024 Waters Corporation. All Rights Reserved.

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