

## Waters Breeze QS HPLC System Software Upgrade Program

### ENHANCE THE CAPABILITIES OF YOUR BREEZE QS HPLC SYSTEM

Enhance the capabilities of your Breeze™ QS HPLC System by taking advantage of the Empower™ QS Software Upgrade program. The upgrade offers greater method flexibility and it keeps your laboratory compatible with new developments in Waters™ HPLC detectors for all your future analytical requirements.

Empower QS Software provides users with an easy-to-use, single-window interface for acquiring, processing, and reporting chromatographic results. It simplifies the complete management of your chromatographic data, providing the information you need at your fingertips, so you can efficiently make informed and accurate decisions.



*The Breeze QS HPLC System continues the balance between Waters world-class technology, high performance, simplicity, reliability, and affordability.*

### THE EMPOWER QS SOFTWARE UPGRADE PROVIDES:

- Microsoft® Windows® 10 operating system support
- Increased productivity and minimal training requirements
- Harmonized control of and streamlined access to a wide variety of analytical chromatographic and detection techniques
- Improved speed and accuracy of laboratory results

### INSTRUMENT PLAN CUSTOMERS - SPECIAL OFFER

- Empower QS Software Support and Maintenance Plan at 50% Discount
- Computer at 20% discount
- Free Installation Services

## UPGRADE OPTIONS

You can choose from two upgrade options to suit your laboratory needs

Software Upgrade Options		
	Empower QS Software Upgrade Only	Empower QS Software, Computer with Installation
<b>Instrument Plan Customers</b> (no SW plan)	<b>50% off SW plan purchase</b> <i>(provides SW upgrade and GPC option plus one-year maintenance/support)</i>	<b>20% off computer list price and free installation</b>
<b>Non-Plan Customers</b>	<b>Purchase SW plan/upgrade</b> <i>(provides SW upgrade and GPC option plus one-year maintenance/support)</i>	<b>20% off computer list price and 50% off installation</b>

### INSTRUMENT PLAN CUSTOMERS:

#### Special Money-Saving Offer

If your system is covered under a current instrument service plan or warranty, you can purchase the Empower QS Software Support and Maintenance Plan at 50% off plan list price. If you require the Breeze QS Software Upgrade with a fully-configured computer, you can purchase separately at a promotional price of 20% off list price. Also, as part of your special offer, Waters will provide upgrade installation services FREE of charge.

### NON-PLAN CUSTOMERS

If your system is not covered under an instrument service plan or warranty, you can purchase the Empower QS Software Support and Maintenance Plan in order to upgrade your system's software. As a special offer, you can purchase the computer at 20% off list price and 50% off installation services.

## INSTALLATION SERVICES

Upon receipt of your Empower QS Software Upgrade package, Waters recommends that you schedule an appointment with one of our certified Field Service Experts to perform the upgrade installation. Doing so ensures all is operational and ready to go. Installation will take approximately 2–3 hours.

## SUPPORTED INSTRUMENTS AND COMPUTER CONFIGURATION SPECIFICATIONS

By upgrading to Empower QS Software, you receive harmonized control of and streamlined access to a wide variety of analytical chromatographic and detection techniques and continuous support for your existing system.

Empower QS is a compliance-ready software allowing you to easily utilize our Qualification Services for all your modules either through our SQT or QWB.

The following instruments are supported:

PUMPS	DETECTORS	INJECTORS	OPTIONAL COMPONENTS
Waters 1525 Binary Pump	Waters 2489 UV/Visible Detector	Manual Injectors	1500 Series Column Heaters
Waters 1515 Isocratic Pump	Waters 2998 Photodiode Array (PDA) Detector	and/or	Extended Temperature Range Column Heater
1525 $\mu$ and 1525EF Pumps*	Waters 2414 Refractive Index (RI) Detector	Waters 2707 Autosampler	eSAT/IN™ to add third-party detectors
	Waters 2475 Multi- $\lambda$ Fluorescence Detector	Waters 2707 Autosampler with sample cooling	1500 Series Vacuum Degassers
	Waters 2424 Evaporative Light Scattering (ELS) Detector		

\* via Ethernet proxy under PCM II.

## Computer hardware configuration available from Waters:

- Lenovo ThinkStation P-520
- Xeon W-2125 4.0 Ghz 4 Cores
- 4 x 16 GB DDR4 2666 ECC RDIMM
- Quadro P1000 4 GB 4 x mini DP (graphics card)
- 2 x 6 TB HD 7200RPM 3.5" SATA
- 1 x 512 GB SSD M.2 PCIe OPAL
- 1 x 10TB HD 7200RPM 3.5" SATA3
- 1 x I210-T1 Single Port Gigabit ET
- 1 x DVD CD-RW

### Operating System

- Microsoft Windows 10

Please refer to the Breeze QS HPLC System and Software Installation and Configuration Guide for more details. Note that Microsoft® Windows® 10 operating system is required.

## INSTRUMENT PLAN FOR BREEZE QS HPLC SYSTEM

Waters Instrument Plans provide you with the tools and support you need to keep your systems performing at the highest standards of performance and reliability.

With a Breeze QS HPLC System Instrument Plan, you get:

- Priority on-site response time
- Priority response telephone support
- Annual Performance Maintenance visit (includes a Performance Maintenance kit and installation)
- Labor and travel for on-site repairs
- Genuine Waters Quality Parts™
- Discount on hardware upgrades
- Discount on customer training
- Firmware updates

## EXCLUSIONS

Normal wear and tear maintenance parts are customer installable. The parts are provided at no additional charge, while the labor and travel are billable events when performed outside of the annual Performance Maintenance on-site visit.

## REQUESTING YOUR UPGRADE IS EASY

To request your upgrade, or for more information, go to [www.waters.com/softwareupgrades](http://www.waters.com/softwareupgrades) or contact your local Waters representative.

Note to Instrument Plan Customers: In order to provide you with the software upgrade, please have your support ID number available to validate your entitlement and ensure appropriate pricing and discount offers.

For your local sales office, please visit [www.waters.com/contact](http://www.waters.com/contact).

# Waters

THE SCIENCE OF WHAT'S POSSIBLE.™

Waters, The Science of Whats Possible, Empower, Breeze, and Waters Quality Parts are trademarks of Waters Corporation. All other trademarks are the property of their respective owners.

©2019 Waters Corporation. Produced in the U.S.A. September 2019 720002868EN LM-PDF

## WATERS GLOBAL SERVICES

### Delivering world renowned services and support

Waters Service and Support offerings are tailored to optimize your laboratory productivity while addressing your budget realities. Our offerings help maintain system peak performance, minimize down time, address scientific application challenges, and support stringent compliance requirements. As your services and support provider, we are committed to the success of your laboratory and business.

Waters quality support and consultative services ensure your success wherever your laboratory is located in the world.



Waters has consecutively earned the ACE Award since 2001 for providing best-in-class technical knowledge, issue resolution, and process support.

Waters Corporation  
34 Maple Street  
Milford, MA 01757 U.S.A.  
T: 1 508 478 2000  
F: 1 508 872 1990  
[www.waters.com](http://www.waters.com)